



Automatic Overbed Table/Model 4528
Acute Care Tables

Model 4528**Top**

Walnut wood grain high pressure laminate on 3/4" particle board core. Flush mounted vinyl T-mold. 18"W x 32"L.

Frame

Chrome plated 11 and 12 gage welded tubular steel column and base.

Height

30"-45" adjustable in 15 one inch increments.

Base

Heavy duty U-base with 28" clearance for easy use with chair/wheelchair. 1 1/2" x 2 1/2" tubular steel. Column mounts on right or left side.

Casters

2" twin wheel



The **AmFab** Deluxe Automatic Lift Overbed Table combines durable, heavy duty construction and attractive functional design with effortless counterbalanced spring height adjustment.

The Deluxe Automatic Lift Overbed Table raises with the slightest upward pressure applied to any point under the top, and locks into 15 positions between 30" and 45" high, and lowers with a touch of the convenient release lever.

The extra heavy duty "U" base has 28" of inside clearance to accommodate most wheelchairs, geri-chairs, and patient room seating. The large 18" x 32" top surface provides plenty of room for eating, reading, or any activity and will support up to 100 pounds evenly distributed.

The optional pullout vanity, Model 8400, features two compartments with hinged acrylic mirrors.

A variety of optional high pressure laminate tops are offered to complement your room decor and remain attractive for years of service. They are designed to resist damage from spills of alcohol, boiling water, food acids, or inks.

Weight

59 lbs. Cube 2.8. UPS Shippable.

Shipping

One table per carton, knocked down. Table easily assembled with four bolts.

Options (Available at additional cost)

- Alternate laminates*
- Vanity & mirror (factory installed when ordered with the table)
- Painted frame*

**Please allow an additional six weeks for delivery.*

AMFAB™



1446 South 35th Street
Galesburg, MI 49053

Phone 269.665.6703
800.432.6322
Fax 269.665.6717