


Bovie DISPOSABLE LOOPS, STERILE

Description: All Bovie LLETZ Loop and Square electrodes feature tungsten wire for superior shape and integrity throughout the excision procedure. They use a standard 3/32" (2.3 mm) shaft. All molded loops are packaged in a patented protective tray to prevent damage during shipping. Ball electrodes are available in 3 mm & 5 mm sizes for fulguration and desiccation while the needle electrodes are available for pinpoint coagulation.



Disposable Loops, Sterile



ES08 5 mm X 5 mm loop – 5/box



ES09 10 mm X 10 mm loop – 5/box



ES10 15 mm X 8 mm loop – 5/box



ES11 15 mm X 10 mm loop – 5/box



ES12 20 mm X 8 mm loop – 5/box



ES13 5 mm X 5 mm loop – 5/box



ES31 20 mm X 20 mm loop – 5/box



ES42 20 mm X 12 mm loop – 5/box



ES43 15 mm X 15 mm loop – 5/box



ES44 15 mm X 5 mm loop – 5/box



ES45 13 mm X 13 mm loop – 5/box



ES46 10 mm X 5 mm loop – 5/box



ES47 25 mm X 10 mm loop – 5/box



ES49 20 mm X 10 mm loop – 5/box



ES51 13 mm X 8 mm loop – 5/box



ES52 10 mm X 8 mm loop – 5/box



ES53 15 mm X 12 mm loop – 5/box

Disposable Squares, Sterile



ES14 5 mm X 5 mm loop – 5/box



ES15 20 mm X 4 mm loop – 5/box




ES16 10 mm X 8 mm loop – 5/box



ES41 5 mm X 10 mm loop – 5/box

Bovie REUSABLE LOOPS, NON-STERILE

Description: All Bovie Reusable LLETZ Loop and Square electrodes are Autoclavable and feature tungsten wire for superior shape and integrity throughout the excision procedure. They use a standard 3/32" (2.3 mm) shaft. All reusable molded loops are individually packaged in a foam insert to prevent damage during shipping. Ball electrodes are available in 3 mm & 5 mm sizes for fulguration and desiccation while the needle electrodes are available for pinpoint coagulation.



Reusable Loops, Non-Sterile



ES08R 5 mm X 5 mm loop – 1/box



ES09R 10 mm X 10 mm loop – 1/box



ES10R 15 mm X 8 mm loop – 1/box



ES11R 15 mm X 10 mm loop – 1/box



ES12R 20 mm X 8 mm loop – 1/box



ES13R 5 mm X 5 mm loop – 1/box



ES31R 20 mm X 20 mm loop – 1/box



ES42R 20 mm X 12 mm loop – 1/box



ES43R 15 mm X 15 mm loop – 1/box



ES44R 15 mm X 5 mm loop – 1/box



ES45R 13 mm X 13 mm loop – 1/box



ES46R 10 mm X 5 mm loop – 1/box



ES47R 25 mm X 10 mm loop – 1/box



ES49R 20 mm X 10 mm loop – 1/box



ES51R 13 mm X 8 mm loop – 1/box



ES52R 10 mm X 8 mm loop – 1/box



ES53R 15 mm X 12 mm loop – 1/box

Reusable Squares, Non-Sterile



ES14R 5 mm X 5 mm loop – 1/box



ES15R 20 mm X 4 mm loop – 1/box



ES16R 10 mm X 8 mm loop – 1/box



ES41R 5 mm X 10 mm loop – 1/box