



Examination light

CoolSpot II

The advanced healthcare lighting solution that fits both your space and your budget.

The aesthetically appealing CoolSpot II offers incredibly pure, white light along with a guaranteed no-drift arm, high-intensity setting and variable spot size control. The CoolSpot II offers intensity of up to 60,000 lux (5574 fc) at 1 meter. The CoolSpot II examination light is the spotlight of choice.

Features

- Up to 60,000 lux (5574 fc) at 1 meter
- 4300 K color temperature
- CRI (Color Rendering Index) of 94
- Bi-level dimmer for variable intensity
- Variable spot size of 2" - 10" (5 cm - 25 cm)
- New, improved and easy to clean housing design with convenient SteriHandle for positioning
- Removable, autoclavable center SteriHandle design allows removal or replacement of SteriHandle with one hand (for ensuring sterile procedures)
- 150-watt halogen bulb provides 600-hour average bulb life on standard setting
- Friction knob allows user to adjust friction on the self-balancing arm positioning to desired tension without a tool
- Mounting systems meet California Seismic Codes
- 120V model standard; 230 / 240V models also available
- UL/IEC 60601-1, IEC 60601-1-2, EN/IEC 60601-2-41 and CAN/CSA C22.2 NO. 601.1 M90 certified
- 5-year limited warranty
- Assembled in USA



Ordering information

120V Model	230/240V Model	230/240V SP Model*	100V Model	Description
CS316FL	CS325FL	C302FL25	CS315FL	Floorstand
CS316SC	CS325SC	C302SC25	CS315SC	Single ceiling mount
CS316DC	CS325DC	C302DC25	CS315DC	Double ceiling mount
CS316W	CS325W	C302W25	CS315W	Wall mount
CS316ST	CS325ST		CS315ST	Single Fastrac mount
CS316DT	CS325DT		CS315DT	Double Fastrac mount
Accessories				
0007006PK				CoolSpot II replacement bulbs, 3/box
1008466B				CoolSpot II handle replacement
0008100PK				Light handle covers, sterile, 25/box
1008259				42" drop tube
1008060				Special length, max. 58" drop tube

*SP Model refers to Spanish Language Manuals.