

# MedWurx



## Stainless Steel Instrument Tables

### Instrument Back Table



#### Features

- Type 304 Solid Stainless Steel Construction with Micro-Satin Finish for Years of Use
- Designed and Manufactured for Durability and Ease of Use in the Operating Room
- 4" Standard, Non Marring Swiveling Easy Roll Casters; 2 Locking
- Optional Upgrade To Antimicrobial Casters (2-Swivel, 2-Locking) , Lips and Drawers
- All Shelves are Hemmed and Sealed for Complete Drip Edge
- Heavy Duty 1.5" EZ Clean Square Stainless Steel Legs for Ease of Cleaning, Strength as well as Durability

#### Product Specifications

- Item ID: MW5115
- 60" L X 24" W X 34" H
- Limited Lifetime Warranty on Stainless Steel
- All Tables Come with Lower Shelf
- Additional Sizes Available, Including
  - MW5112 20" L X 16" W X 34" H
  - MW5113 36" L X 24" W X 34" H
  - MW5114 48" L X 24" W X 34" H
  - MW5115 60" L X 24" W X 34" H
  - MW5116 72" L X 24" W X 34" H
  - MW5313 36" L X 20" W X 34" H
  - MW5314 48" L X 20" W X 34" H
- Many Options and Models Available - see Website for Latest Details
- Made in the USA

rev311220

**www.MedWurx.com**