




### Advantages

- ✓ **EMR Connectivity**  
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- ✓ **Motion Detection**  
Scale has the ability to weigh accurately despite movement of the subject on the scale
- ✓ **Value-Added Features**  
Operating functions (LB/KG Conversion, KG Lock Out, Zero Out/Tare, Hold/Release, Auto Off) included at no additional cost
- ✓ **Measuring Tape Included**  
Easy-to-read measuring tape allows user to measure both weight and height simultaneously



| Model | Shipping Dims. (L x W x H) / Cube | Weight (LB) | UPC Code     |
|-------|-----------------------------------|-------------|--------------|
| 553KL | 28.125" x 19" x 5.5" / 1.70       | 12.25 lb    | 892076002329 |



### Specifications

- Capacity: 44 lb / 20 kg
- Resolution: 0.5 oz / 10 g
- Tray Size: 24 7/8" (w) x 13" (d) x 1 3/8" (h) / 632 mm (w) x 330 mm (d) x 35 mm (h)
- Measuring Tape: 0" - 24" / 0 cm - 58 cm
- Display: 1" LCD
- Power Source: 6 AA batteries included, optional power adapter (order # ADPT30)
- Connectivity: USB
- Carrying Case: Optional carrying case (order # 553CASE)
- Functions: LB / KG Conversion, LB / KG Lock Out, Zero Out / Tare, Hold / Release, Auto Off
- 2 Year Limited Warranty



Carrying Case for 553KL  
Item # 553CASE