

Nicolet *Elite*TM - OB applications

The *Elite* incorporates all of the best features of over 20 years of Imex Doppler experience. Mix and match these features to create your Doppler of choice!



- Optional digital heart rate display to easily determine fetal heart rate.
- Standard 9V alkaline or optional rechargeable batteries available.
- The comfortable shape of the *Elite* probe is easy to hold and features broad beam technology, making it very easy to find and lock onto the signal of interest.
- The new probes also incorporate Static Suppression Doppler (SSD) technology for a clear, static-free signal and the best possible diagnostic information.
- New! 2MHz waterproof probe is available for underwater labor and delivery. The main unit remains outside of the water bath to provide clear, undistorted sound.
- Obstetrical applications include detecting pregnancy, reassuring the mother, and determining multiple pregnancies.

For more information on *Elite* model configurations and to compare to existing Nicolet products visit www.medicalresources.com or call 1-800-860-4716

Nicolet *Elite* Technical Specifications

Dimensions:	6.8" x 2.9" x 1.0" (17cm x 7.5cm x 2.5cm)
Weight:	0.60 lb. (0.25kg)
Power source:	110V recharger (220 V optional): Models E100R, E200R; Battery powered: Models E100, E200
Battery type:	8.4 Volt NiCad user replaceable: Models E100R, E200R; 9V alkaline: Models E100, E200
Battery recharge time:	10 hours (3.5 hours to restore 50% of capacity): Models E100R, E200R
Display range:	60-199 +/- 2: Models E200, E200R
Display Resolution:	1 beat per minute: Models E200, E200R