



647
646

MIDMARK 647 AND 646 PODIATRY CHAIRS

Rediscover the most important part of the podiatry chair, the foot section.



Ergonomic Foot Section
The easy-to-reach, release handles allow for one-handed, smooth, quiet operation. Integrated touch controls means there is no hand control to get knocked to the ground, eliminating costly breakage and allowing unrestricted patient access.
Patent pending



Chair Arm Rotation
Chair arms offer the patient a sturdy platform for accessing the chair. The arms rotate completely out of the way, providing accessibility for patients with mobility concerns, such as disabled or obese patients.



Aseptic Design
Covered glides on the foot section and debris tray help keep dirt and debris out to protect the smooth operation and to provide a clean, aseptic working environment your patients will appreciate.
Patent pending



Chair Arms
Chair arms can be easily installed or removed for your convenience...all without the use of tools.



Sealed Foot Control
Midmark's foot control is sealed against dirt and debris, making cleaning simple and ensuring long life. Quick Chair™ button returns the table to the chair position for quick patient exit.



Quick Exam™
The Midmark 647 non-programmable foot control offers a Quick Exam button that moves the tilt and seat simultaneously, repositioning the chair to the exam position in one smooth motion.

Powered Multi-Positioning
The Midmark 647 offers a Barrier-Free® low height of 19", allowing disabled and obese patients easy accessibility. Our custom designed motors quickly move your patients into proper exam position, allowing you to be more productive. Powered movement of the seat, back, and tilt gives you the flexibility of positioning patients where you need them.



low height position



flat position



exam position

The Midmark 647 Barrier-Free® Podiatry Power Procedures Chair offers more standard features.

The 647 podiatry chair has a host of options that add functionality, versatility and improve efficiency.



With multiple options available, you can purchase a chair that is custom-suited for you and your practice, delivering unsurpassed efficiency and functionality for a wise, long-term investment.



Rotation
Chair rotation can provide you with the needed space to work, even in the smallest of rooms. A simple touch of the pedal on either side of the chair permits 350° of rotation; touch the pedal again to lock the chair in place. It's that easy.



Programmability
Optional One-Touch, programmable positioning quickly moves all sections of the chair to the programmed position simultaneously...speed, efficiency and simplicity with the touch of a button.



Electrical Receptacles
Optional hospital grade electrical receptacles allow easy access to electrical power, freeing your floor space from power cords, creating a clean look in the procedure room.

Accessories that add efficiency and functionality



Disposable Touch Control Covers
Disposable touch control covers make cleanup easy and allow you to maintain an aseptic environment during the procedure. Simply press them on, then peel them off.



Casters
Locking casters provide mobility without sacrificing valuable floor space. Casters add 3 ½" to the minimum and maximum seat height, creating space under the chair for Mayo stand legs.



Vision Block Screen
The vision block screen is simple to install, remove or rotate, allowing you to easily block the patient's view while you perform sensitive procedures. The screen material is color matched to the upholstery of the chair, enhancing aesthetics.



Double Arm Instrument Tray
Place instruments and supplies where needed most...close to the patient's feet. The instrument tray attaches to either side of the chair, saving valuable floor space. Double arm design allows for a wider range of motion. (Single arm version is available.)



Wireless Foot Control
Remove cords and eliminate obstacles. Enjoy the freedom of locating the foot control anywhere in the room. Low battery audible indicator provides early warning of approaching battery replacement.



Base Rails
Base rails provide a convenient mounting point for attaching other accessories. Tucked under the seat, they won't interfere with patient accessibility.



Paper Roll Holder
The ability to cover the back and seat section with paper can be critical if you want your chair to also work for your non-podiatrist partners. Holds 21" x 3 ½" paper roll.
**only available on 647*



Midmark 646 Basic Podiatry Procedures Chair.

The Midmark 646 Basic Podiatry Procedures Chair offers an economical solution to the podiatrist who works exclusively from the seated position. It offers powered back and tilt with a fixed seat height of 21". Touch controls and sealed foot control are standard on this model, as are covered glides on the foot section and debris tray. You'll be surprised at how easy it is to keep your chair clean and operating smoothly.

Powered tilt moves the foot section into the exam position and the foot section extends to accommodate taller patients. The Midmark 646 can also be customized with any of the available accessories, with the exception of the paper roll holder.

647 Specifications

Patient Load Rating:
450 lb. (204 kg)

Base Height:
Minimum 19" (48.3 cm)
Maximum 40" (101.6 cm)

Foot Section Height:
With foot section retracted
15.5" to 54" (39.4 cm to 137.2 cm)

With foot section fully extended
14" to 57" (35.6 cm to 144.8 cm)

Foot section:
26" W x 6" H (66.0 cm x 15.2 cm)

Foot extension: 9" (22.9 cm)

Debris tray: 15.1" W x 10.6" H x 1.5" D
(38.4 cm x 26.9 cm x 3.8 cm)

Back Support Angle Range:
0° to 80°

Seat Tilt Angle Range: 0° to 30°

Seat Width with Chair Arms Installed:
37.5" (95.3 cm)

Foot Control:
14.5" W x 6" H (36.8 cm x 15.2 cm)

Optional Electrical Outlet Rating:
115 VAC, 3 Amps

Options

Premium Upholstery Top:

Seamless design for easy cleaning

Ultra-Premium Upholstery Top:

Stitched, plush design for high-end appearance and unparalleled comfort

Programmability:

Three positions that you customize to your needs

Electrical Outlets:

Duplex, hospital grade
(on each side)¹

Rotation:

Allows 350° of table rotation and adds 3.5" to the minimum and maximum seat height specifications

646 Specifications

Patient Load Rating:
450 lb. (204 kg)

Base Height:
Fixed height of 21" (53.3 cm)

Foot Section Height:
With foot section retracted
18" to 35" (45.7 cm to 88.9 cm)

With foot section fully extended
15.5" to 37.5" (39.4 cm to 95.3 cm)

Foot section:
26" W x 6" H (66.0 cm x 15.2 cm)

Foot extension: 9" (22.9 cm)

Debris tray: 15.1" W x 10.6" H x 1.5" D
(38.4 cm x 26.9 cm x 3.8 cm)

Back Support Angle Range:
0° to 80°

Seat Tilt Angle Range: 0° to 30°

Seat Width with Chair Arms Installed:
37.5" (95.3 cm)

Foot Control:
14.5" W x 6" H (36.8 cm x 15.2 cm)

Optional Electrical Outlet Rating:
115 VAC, 3 Amps

Accessories

Base rail

Caster base¹

Casters¹

Cord, power, 3' length

Disposable touch control covers⁴

Instrument tray, single arm³

Instrument tray, double arm³

IV pole³

Paper roll holder

Vision block screen

Wireless control, linear-style²



Premium Comfort System™

Advanced elastomeric material stretched over a seat cut-out enhances patient comfort for a “chair-like” feel. The wide seat and back provides your patients with the comfort and support they need.

Length measures 75" ...
can be increased to 84"
by extending foot section.

Patent No. US 7,376,991 B2

Ultra-Premium upholstery shown.

¹ – Not available on rotation models

² – Includes foot control
(with batteries)

³ – Requires base rail accessory

⁴ – Roll of 600 sheets

Premium Solid Upholstery Colors

Blueberry

Cashmere Blue

Silver Sage

Terra Cotta

Fossil Grey

Tea Green

Black

Premium Pattern Upholstery Colors

Belagio

Capri

Siena

Treviso

Firenze

Tuscany

Ultra-Premium Solid Upholstery Colors

Ultra-Pacific Blue

Ultra-Sky

Ultra-Spruce

Ultra-Canyon

Ultra-Arctic

Ultra-Olivine

Ultra-Raven Wing



Midmark is an ISO 13485 and ISO 9001 Certified Company
For more information or a demonstration, contact your Midmark dealer or Call: 1-800-MIDMARK Fax: 1-800-365-8631
Outside the U.S.A. Call: 1-937-526-3662 Fax: 1-937-526-8392
or visit our website at midmark.com

© 2008 Midmark Corporation
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286
Products subject to improvement changes without notice
Litho in U.S.A. 007-0730-00 Rev. A1 (9/08)

