

Stainless Steel Blanket Warmer

CounterTop - Solid Door



Features

- Type 300 Solid Stainless Steel Double Walled and Insulated Construction
- Antimicrobial Handle
- 4.5 Total Cubic Feet
- One Adjustable Perforated Stainless Steel Shelf
- Self-Closing Energy Efficient Door with Non-Latex Magnetic Gaskets
- Forced Air Circulation for Consistent Ambient Air Temperatures
- Illuminated On \ Off Switch
- Programmable and Visual Alarm Digital LED Controls (Operates in Both F and C)

Product Specifications

- Item ID: MW01121
- Exterior Dimensions: 30" W X 24" D X 24" H
- Interior Dimensions: 28" W X 17.5" D X 16" H
- Electrical is 120V (optional 220V available), 50/60 Cycles
- Blanket Warmers have a Limited Lifetime Warranty on Stainless Steel, Three Year Limited Warranty on All Parts and One Year Limited Warranty on Labor
- Models and Options Available: Single and Double Compartment, Full Size, Under Counter, Counter Top, Solid Doors, Glass Doors, Countertop Stand and Mobile Base
- Many Options and Models Available see Website for Latest Details
- ETL Listed with Hospital Grade Plug and Limit Switch
- Made in the USA rev317011

www.MedWurx.com