

MedWurx

Stainless Steel Blanket Warmer

Full Size - Double Glass Doors



Features

- Type 300 Solid Stainless Steel - Double Walled and Insulated Construction
- Antimicrobial Handles
- 15.1 Total Cubic Feet
- Two Adjustable Perforated Stainless Steel Shelves per Compartment
- Self-Closing Energy Efficient Doors with Non-Latex Magnetic Gaskets
- Forced Air Circulation for Consistent Ambient Air Temperatures
- Illuminated On \ Off Switch
- Programmable and Visual Alarm Digital LED Controls (Operates in Both F and C)

Product Specifications

- Item ID: MW01222
- Exterior Dimensions: 30" W X 24" D X 72" H
- Interior Dimensions: Upper : 28" W X 17.5" D X 24" H
Lower: 28" W X 17.5" D X 30" H
- Electrical is 120V (optional 220V available), 50/60 Cycles
- Blanket Warmers have a Limited Lifetime Warranty on Stainless Steel, Three Year Limited Warranty on All Parts and One Year Limited Warranty on Labor
- Many Options and Models Available - see Website for Latest Details
- Models and Options Available: Single Compartment, Full Size, Under Counter, Counter Top, Solid Doors, Glass Doors, Countertop Stand and Mobile Base
- ETL Listed with Hospital Grade Plug and Limit Switch
- Made in the USA

rev317011

www.MedWurx.com