

Ritter 230 Specifications

Patient Load Rating: 450 lb. (204 kg)
Height: Minimum 22.5" (57.2 cm)
Maximum 40" (101.6 cm)
Seat Tilt Angle Range: 0° to 30°
Back Support Angle Range: 0° to 85°
Foot Support Angle Range: 0° to 90°
Upholstery Width Dimensions: Headrest: 18" (45.7 cm)
Back support: 27" (68.6 cm)
Seat: 27" (68.6 cm)
Foot support: 18" (45.7 cm)
Patient Support Surface Length:
Flat position: 72" (182.9 cm)
Flat position with headrest extended: 83" (210.8)
Paper Roll Holder Capacity: Holds 18" x 3.5" paper roll (45.7 x 8.9 cm)
Stainless Steel Treatment Pan Dimensions: 11.5" L x 9" W x 2.5" D
(29.2 x 22.9 x 6.4 cm)
Optional Electrical Outlet Rating: 115 VAC, 3 Amps

Options

Electrical Outlets: Duplex, hospital grade (on each side)