



Executive Split-top Overbed Table / Model 4850H, 4850U

## Acute Care Tables



Executive Split-top Overbed Table / Model 4850H, 4850U

Model 4850U

**Slide activation tabs**

2 slide activation tabs are located on column side of the tabletop. The location of the tabs provides patient and caregiver convenience while also promoting safety.

**Height**

Automatic gas cylinder lift mechanism adjusts height smoothly and quietly.

**Frame**

Base & columns are powder coated 11 and 12 gauge welded tubular steel.

**Top**

When the split-top table is activated, the 18" x 32" upper table slides either forward or back allowing the 17" x 26" lower table to move into the opposite position. This allows for more user-friendly work or eating space. The high-pressure laminate top, with resilient T-molding edge, resists chipping, scratching and stains.

**Base**

Constructed of 1 1/2" x 2 1/2" tubular steel. Available in Heavy-duty U (shown) or H-base. Bases have 28" clearance to provide proper space for the positioning of wheelchairs and other similar seating systems. Column can be mounted on either right or left side.

**Casters**

2" twin wheel casters provide smooth maneuvering.

Optional factory installed vanity features 2 way de with non breakable acrylic mirrors and rounded corners for easy cleaning.

The **AmFab** split-top table features the reliable gas cylinder lift mechanism for smooth and quiet height adjustment.

The top automatically rises with the slightest pressure on any point of the table, lowers with the slightest touch of the trip handle.

Available in Heavy-duty U (shown above) or H-base design. The bases have 28" clearance to provide proper space for the positioning of wheelchairs and other similar seating systems. The generous 18" X 32" top surface provides plenty of space for eating, reading or light therapy and will support 100lbs evenly distributed. The user-friendly lower 17" x 26" table is an excellent extra space for food trays or activities while the upper top is in use.

The optional pullout vanity, 8400, provides space for personal items and features two compartments with hinged acrylic mirrors.

A variety of optional high-pressure laminate tops are offered to compliment any room décor and will remain attractive for years of service. They are designed to resist damage from spills of alcohol, boiling water, food acids or inks.

**Weight** 75 lbs., Cube 12.8

**Shipping**

One table per carton, assembled, ships via Common Carrier.

**Options** (Available at additional cost)

- Alternate laminates
- Vanity & Mirror (factory installed)
- Alternate frame colors

**AMFAB™**



1446 South 35th Street  
Galesburg, MI 49053

Phone 269.665.6703  
800.432.6322  
Fax 269.665.6717