

## Mobile Adult Mirror with Cuff Weight Rac

### MODEL 6210-9

- 1/4" thick distortion-free glass mirror with ANSI 297 safety backing
- All laminate construction
- Rolls easily on 3" swivel casters
- ABS plastic bumpers on base sides
- Cuff weight sliding hook system makes lateral positioning fast and easy. 12 cuff weight hooks included
- Holds up to 48 cuff weights (Not included)



California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals including 1, 3-Butadiene which is known to the State of California to cause cancer and reproductive harm. For more information go to: [www.P65Warning.ca.gov](http://www.P65Warning.ca.gov).

7



### MODEL 6210-9 SPECIFICATIONS

|         |                              |
|---------|------------------------------|
| Overall | 31.5" W x 20.5" D x 75.75" H |
| Base    | 31.5" W x 20.5" D            |
| Top     | 28.25" W x 4" D x 70.75" H   |

### Laminate



Maple