

## Mobile Adult Mirror

### MODEL 6210


- 1/4" thick distortion-free glass mirror with ANSI 297 safety backing
- All laminate construction
- Rolls easily on 3" swivel casters
- ABS plastic bumpers on base sides

#### Shipping Information:

- Ships on a pallet that is 139 pounds and 80" x 38" x 18"
- The mirror alone weighs approximately 85 pounds



California residents please be advised of the following, pursuant to Proposition 65:

**7**  
 **WARNING:** This product can expose you to chemicals including 1, 3-Butadiene which is known to the State of California to cause cancer and reproductive harm. For more information go to: [www.P65Warning.ca.gov](http://www.P65Warning.ca.gov).



#### MODEL 6210 SPECIFICATIONS

Overall	31.5" W x 20.5" D x 75.75" H
Base	31.5" W x 20.5" D
Top	28.25" W x 4" D x 70.75" H

#### Laminate



Maple