

## Fun Series, Frosty Friends Ice Cream Truck, Treatment Table



FRONT & DOOR SIDE



BACK & OPPOSITE SIDE

### MODEL 7070

- Clinton original illustrations
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Graphics are on two long side and two short ends, manufactured as pictured
- Optional locking doors available
- All colors are as shown
- 3 Storage compartments with one adjustable shelf each
- Soft-close, door hinges
- Chip resistant PVC edges
- Seamless rounded corner top
- Premium, stain resistant, woven knit back upholstery (black)
- Adjustable feet
- Doors are only on one side
- Paper dispenser included
- Special order reverse side doors are available
- 350 lbs. load capacity under normal use

### 7070 SPECIFICATIONS

72"L x 27"W x 34"H

California residents please be advised of the following, pursuant to Proposition 65:

**!** **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).


1

MANUFACTURED IN THE  
**USA**

### Options




066 Door & Inside Latch Combo




055 Door Locks



040 Paper Cutter



20 Small Pillow



60 Large Pillow

### Upholstery



Black