


## Multi-Use, Ultrasound Table with Stirrups



### MODEL 80069

- Three-piece top with pneumatic adjustable backrest, manual footrest, and stirrups
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Trendelenburg/Reverse Trendelenburg plus or minus 15 °
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Industry leading 5-year, limited warranty
- Optional side rails & casters available
- 34" wide model available
- Model 80069: 500 lbs load capacity under normal use
- Model 80069-X: 600 lbs load capacity under normal use

California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).



*Trendelenburg*



*Reverse Trendelenburg*



*Optional Side Rails & Casters*

### 80069 SPECIFICATIONS

	length	width	height range
overall	76"	27"	21 <sup>1</sup> / <sub>2</sub> " – 37 <sup>1</sup> / <sub>2</sub> "
backrest	38"	27"	
footrest	12"	18 <sup>1</sup> / <sub>2</sub> "	

80069-X Extra-Wide 34" width model



### Upholstery Colors



## OPTIONS

### Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set

080	length	width
	92"	2"
	233.67 cm	5 cm



### Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

098



### Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087



### Paper Holder

- Clear plastic with elastic strap and buckle

040

