

Flat Top, Imaging Table with Drop Window



MODEL 80071

- One piece flat top with drop window
- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Industry leading 5-year, limited warranty
- Model 80071: 500 lbs load capacity under normal use
- Model 80071-X: 600 lbs load capacity under normal use

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

80071 SPECIFICATIONS

	length	width	height range
overall	72"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "
window	10"	10"	

80071-X Extra-Wide 34" width model



Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

OPTIONS

Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set

080	length	width
	92"	2"
	233.67 cm	5 cm

Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

098

Power Backrest

- Smooth, power control raises to 80° and lowers to flat

011

Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087

Paper Holder

- Clear plastic with elastic strap and buckle

040

