

## Open Base Power Casting Table with ClintonClean™ Leg Rest

### MODEL 80350

- Pneumatic backrest is adjustable from either side and raises 80 degrees from flat
- ClintonClean™, Solid Plastic Surface, leg rest is adjustable from either side and lowers up to 80 degrees from flat
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Seamless, rounded corner, upholstered top
- 2" ultra-firm padding
- Low height may help achieve ADA tax credit
- Customize with options
- Made in the USA by Clinton Industries
- 500 lbs. (227 kg) load capacity under normal use

# POWER 500



### 80350 SPECIFICATIONS

Overall	76" L x 27" W x 18–35" H
Backrest	38" L x 27" W
Footrest	13" L x 26" W

California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1



### Upholstery Colors



Slate Blue

Wedgwood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black

## OPTIONS

### Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set (92" x 2")

080



### Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

082



### Hand-Held Controller

- Alternative to foot control
- Stores easily under top when not in use

084



### Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087



### 030 Paper Dispenser

- Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide


030



### Paper Holder

- Clear plastic with elastic strap and buckle

040



### Paper Dispenser/Holder Combo

- Paper dispenser and paper holder ordered as one option

034