

Open Base, Extra Wide, Bariatric, Power Table with Adjustable Backrest

MODEL 84200

- Extra wide, bariatric model with choice of two widths
- Dual-cylinder, pneumatic backrest is adjustable from either side and raises 80 degrees from flat
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Rounded corner, upholstered top
- 2" (5 cm) ultra-firm padding
- Low height may help achieve ADA tax credit
- Customize with options
- Made in the USA by Clinton Industries
- 600 lbs. (272 kg) load capacity under normal use

POWER 600 BARIATRIC




84200-34 SPECIFICATIONS

Overall 72" L x 34" W x 18-35" H

Backrest 38" L x 34" W

84200-40 Extra-Wide 40" width model

California residents please be advised of the following, pursuant to Proposition 65:

1
 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.



Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

OPTIONS

Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set (92" x 2")

080



Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

082



Hand-Held Controller

- Alternative to foot control
- Stores easily under top when not in use

084



Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087



030 Paper Dispenser

- Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide


030



Paper Holder

- Clear plastic with elastic strap and buckle

040



Paper Dispenser/Holder Combo

- Paper dispenser and paper holder ordered as one option

034