

ENT Chair with Foot Operated Pump

Model 8612


The UMF Medical 8612 offers a manually adjustable backrest, footrest, and foot operated pump for hi-lo movement. Offers a 360 degrees rotation, hinged arms, and 250lbs weight capacity.



Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA



River Rock



Sandy Grey



Mocha



Steel Blue



Midnight Blue



Onyx

Key features at a glance

- › 250lb (113.4kg) patient weight capacity
- › Ultra comfort seamless top with antimicrobial guard
- › Sturdy steel frame
- › Available in all UMF Medical upholstery options
- › 360° rotation with swivel lock
- › 10° elevation from horizontal position
- › Manual foot pedal height adjustment
- › Manually adjustable backrest
- › Manually adjustable footrest
- › Hinged arms
- › 3-year limited warranty

Specifications

Min. Seat Height: 23" (58.4cm)
 Max. Seat Height: 38" (96.5cm)
 Overall Width: 31.5" (80cm)
 Seat Width: 18" (45.72cm)
 Arm Length: 23" (58.4cm)
 Length: 71" (180.3cm) fully extended
 Net Weight: 166lb (75.3kg)
 Shipping Weight: 225lb (102kg)
 Shipping Dimensions: L44xW34xH64 (111.8cm x 86.4cm x 162.6cm)

Accessories

031 – Optional backrest with articulating headrest