

E-Series, Power Table with Adjustable Backrest and Drop Section

MODEL 86300

- Features an adjustable, pneumatic backrest and friction lock drop section which are accessible from either side
- Solid, standard drop section (may be ordered with a nose hole as an option) (Please specify when ordering)
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Seamless, rounded corner, upholstered top
- 1.5" (4 cm) ultra-firm padding
- Low height may help achieve ADA tax credit
- Customize with options
- Made in the USA by Clinton Industries
- 450 lbs. (204 kg) load capacity under normal use

E-Series
POWER TABLES



SPECIFICATIONS

Overall 72" L x 27" W x 18.5 –35.5" H

Backrest 38" L x 27" W


Drop Section 13" L x 27" W

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.



Options



009 Nose Hole — Nose hole in drop section of 86220 or 86300

030 Paper Dispenser — Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

040 Paper Holder — Clear plastic with elastic strap and buckle

034 Paper Dispenser/Holder Combo — Paper dispenser and paper holder ordered as one option

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black