




BOVIE[®] BANTAM|PRO[™]

Electrosurgical Generator



The Bovie Bantam|PRO, A952, is a High Frequency Desiccator with cutting ability. You might call it a split personality. One half of the unit is a high frequency desiccator, an upgrade to our very popular Aaron[®] 940. The other half is a Bovie Generator, giving you 50 watts of cutting power and coag capabilities – plenty of power for most cutting procedure performed in the doctor's office.

High Frequency Desiccator with 50 watts of Cutting capabilities

Two machines in one convenient package

Adjustable in 1/10th watt increments in the first 10 watts, for the most accurate output available

No grounding pad required when working in the desiccation mode



BOVIE® BANTAM|PRO™ A952 ELECTROSURGICAL GENERATOR



STANDARD EQUIPMENT

A952 (1) electrosurgical generator,
A902 (1) handpiece
A910 (2) non-sterile disposable handpiece sheaths*
A910ST (2) sterile disposable handpiece sheaths*
ES01 (3), **ES02** (1), **ES20** (1) sterile disposable electrodes
A804 (5), **A806** (5), **A806DE** (2) disposable non-sterile dermal tips*
A805 (2), **A807** (2), **A807DE** (2) disposable sterile dermal tips*
A1252C (2) reusable grounding cord
ESRE (5) disposable split grounding pads
A837 wall mount kit*
 110 VAC Hospital-grade power cord*
 User's guide (MC-55-238-001)*



A902



ES01, ES02, ES20




A1252C



ESRE

* Not Pictured

COMPLEMENTARY PRODUCTS†



A813
Table Top Stand †



A803
Footswitch †




A827V
Forcep Cable †



A824
Reusable Bipolar Forceps †



Lletz Loop
Electrodes †



Resistick II
Coated Electrodes †



Tungsten Wire
Loop Electrodes †

† Sold Separately

PRODUCT SPECIFICATIONS