

PowerSpin™ DX Centrifuge

Product Images



Short Description

1. POWERSPIN™ DX CENTRIFUGE
2. VARIABLE SPEED 1000-3400 RPM
3. VARIABLE SPEED 500-3500 RPM
4. PRE-PROGRAMMED OPTIONS FOR BLOOD, URINE AND FECAL
5. 6 PLACE HORIZONTAL ROTOR
6. 60 MIN DIGITAL TIMER
7. CAPACITY 6X15 ML TUBES
8. POWER INPUT 110V

Description

UNICO is proud to present the PowerSpin™ DX Centrifuge

With 3 rotor options this centrifuge is a great answer for every day laboratory tasks. Choose from 24 place fixed rotor for 2-15ml tubes (#C8724), 6 place horizontal rotor for 2-15ml tubes(#C8706) or 4 place horizontal rotor for up to 50ml tubes (#C8704)

The PowerSpin DX series features a large, spacious stainless steel lined interior, all digital push-button control with 3 preset operational modes for Blood, Urine and Fecal, variable speed 500-4000 RPM and zero RPM lid latching system. DX's are available in three main versions, fixed angle and horizontal spin rotors. All DX's feature maintenance-free brushless motors and suction cup rubber feet.

Additional Information

Accessories	Header	Accessories
	C8600-01	6 Place Horizontal Rotor, 2-15ml for DX Centrifuge
	C8600-02	4 Place Horizontal rotor, 50ml DX Centrifuge
	C8600-03M	24 Place Fixed Angle Rotor, 2-15ml for DX Centrifuge, Metal (without Tube Shields)
	C800-02	Regular tube shield, pack of 2
	C800-03	Adapter for short tubes, pack of 10
	C800-04	Pediatric or microtube adapter, pack of 2
	C800-05	Tube Lift Adapters, Pk/12
Specification	Speed (RPM)	500~4000 rpm
	Maximum RCF (g)	2508
	Rotor	6 place horizontal
	Tube Size	2-15ml
	Motor	DC
	Power Input	110V
Size and Weight	Size	13.75 inch Height (22 inch height with lid open) x 15 inch Deep x 11 inch Width
	Shipping Size	19"X18"X15"
	Weight	31 LBS
	Shipping Weight	36 LBS