

SIMPLIFY YOUR DAILY WORK
No other system can give you as much in a single platform

Diagnostic Station DS20

**American Technology &
Swiss Precision**



Items	DS20 Version Masimo 3-lead ECG 0A.800000	DS20 Version Masimo 3-lead ECG, Temperature 0A.800005	DS20 Version Masimo 3 & 12-lead ECG 0A.800010	DS20 Version Masimo 3 & 12-lead ECG, Temperature 0A.800015
Display	18.5" Touch Screen	18.5" Touch Screen	18.5" Touch Screen	18.5" Touch Screen
Vital Signs	<ul style="list-style-type: none"> - ECG Wave Form (3-lead) - SpO2 - SpCO (License req.) - Respiration - NIPB (Auto/Manual) - Weight/Height/BMI through External Digital Scales 	<ul style="list-style-type: none"> - ECG Wave Form (3-lead) - SpO2 - SpCO (License req.) - Respiration - Temperature - NIPB (Auto/Manual) - Weight/Height/BMI through External Digital Scales 	<ul style="list-style-type: none"> - ECG Wave Form (3 & 12-lead) - SpO2 - SpCO (License req.) - Respiration - NIPB (Auto/Manual) - Weight/Height/BMI through External Digital Scales 	<ul style="list-style-type: none"> - ECG Wave Form (3 & 12-lead) - SpO2 - SpCO (License req.) - Respiration - Temperature - NIPB (Auto/Manual) - Weight/Height/BMI through External Digital Scales
Resting ECG	No	No	Yes	Yes
Resting Rhythm ECG	Yes	Yes	Yes	Yes
ECG Channels	3	3	3 & 12	3 & 12
Spirometry	Yes (optional accessory)	Yes (optional accessory)	Yes (optional accessory)	Yes (optional accessory)
Database	Yes	Yes	Yes	Yes
Worklist	Yes	Yes	Yes	Yes
Memory	72 Hours of Vital Signs Trends	72 Hours of Vital Signs Trends	72 Hours of Vital Signs Trends	72 Hours of Vital Signs Trends
ECG Interpretation	No	No	Yes	Yes
Spirometry Interpretation (FVC)	Yes	Yes	Yes	Yes
Printer	External Thermal or Regular Printer Available	External Thermal or Regular Printer Available	External Thermal or Regular Printer Available	External Thermal or Regular Printer Available
Digital Scales	<ul style="list-style-type: none"> - Health o Meter, Model 500KL - Detecto, Model APEX Serial - Fairbanks, Model 28142 	<ul style="list-style-type: none"> - Health o Meter, Model 500KL - Detecto, Model APEX Serial - Fairbanks, Model 28142 	<ul style="list-style-type: none"> - Health o Meter, Model 500KL - Detecto, Model APEX Serial - Fairbanks, Model 28142 	<ul style="list-style-type: none"> - Health o Meter, Model 500KL - Detecto, Model APEX Serial - Fairbanks, Model 28142
Mounts	Mobile Carts / Wall Mounts / Desk Stands (Optional)	Mobile Carts / Wall Mounts / Desk Stands (Optional)	Mobile Carts / Wall Mounts / Desk Stands (Optional)	Mobile Carts / Wall Mounts / Desk Stands (Optional)
Connectivity	<ul style="list-style-type: none"> - LAN Connection - Wi-Fi Connection - Extensive EMR compatibility Connection 	<ul style="list-style-type: none"> - LAN Connection - Wi-Fi Connection - Extensive EMR compatibility Connection 	<ul style="list-style-type: none"> - LAN Connection - Wi-Fi Connection - Extensive EMR compatibility Connection 	<ul style="list-style-type: none"> - LAN Connection - Wi-Fi Connection - Extensive EMR compatibility Connection
PC Storage	Schiller Sema Software	Schiller Sema Software	Schiller Sema Software	Schiller Sema Software
Warranty	2 Years	2 Years	2 Years	2 Years



Diagnostic Station DS20

American Technology & Swiss Precision

Optional Sensors and Accessories



3-lead ECG Cable
(2.400211)



Disposable ECG Electrodes (2.155025)



NIBP Adult Cuff
(2.120019)



NIBP Hoses
(2.100746)



12-lead ECG Cable
(2.400071)



Snap Clip Adapter
(2.155032)



Biotabs disposable electrodes (2.155031)



Barcode Scanner
(SR20001R-UR)



SP-250 Spiro-sensor
(2.100022)



Disposable Mouth Pieces
(2.100077)



Calibration Syringe
(2.100037)



Exergen Thermometer
(2.101050)



SpO2 Masimo Sensor
(2.100598)



Masimo Cable
(2.310319)



Trolley for DS20
(2.101066)



Kangaroo Trolley
(2.100642)



Desk Stand
(2.100491)



Wall Mount
(2.100492)



Health o Meter Scale
(500KL)



Detecto Scale
(APEX)



Fairbanks Scale
(28142)



SpCO Masimo License
(4.150290)