


**International:** For complete occupational screening using accurate and easily understood targets to help determine effective vision. Includes Purdue Job Standards Booklet to easily identify individuals requiring further testing for all occupations. Intermediate lenses available for VDT operators, which is the fastest growing population in the industry today. Now test Distance, Intermediate, and Near Point Vision.




Slide #1 (2000-255) Far Point  
Acuity Both Eyes Binocular  
(20/200-20/17)




Slide #7 (2000-254) Far Point  
Lateral Phoria (1 Diopter Increments)




Slide #2 (2000-256) Far Point  
Acuity Right Eye Monocular  
(20/200-20/17)



Slide #8 (2000-260) Near Point  
Acuity Both Eyes Binocular  
(20/200-20/17)




Slide #3 (2000-257) Far Point  
Acuity Left Eye Monocular  
(20/200-20/17)



Slide #9 (2000-261) Near Point  
Acuity Right Eye Monocular  
(20/200-20/17)




Slide #4 (2000-258) Far Point  
Stereo Depth Perception  
(362-9.7 Seconds of Arc)



Slide #10 (2000-262) Near Point  
Acuity Left Eye Monocular  
(20/200-20/17)



Slide #5 (2000-259) Far Point  
Color Perception (Pseudo Ishihara)



Slide #11 (2000-264) Near Point  
Lateral Phoria (1 Diopter Increments)



Slide #6 (2000-253) Far Point  
Vertical Phoria (1/2 Diopter Increments)



Slide #12 (2000-263) Near Point  
Vertical Phoria (1/2 Diopter Increments)