

## **ABS ABT-FFP-04 Details and Specifications**

4 Cu. Ft. capacity

Two fixed shelves with three removable baskets

Keyed door lock

Magnetic door gasket for positive seal

Integrated door handle

HFC-free refrigerant (R600a), foam insulation and packaging

Two year parts and labor warranty

*Note: This unit must have 4" clearance on sides and back for adequate ventilation*

Units meet NFPA guidelines 45 and 70

Microprocessor temperature controller

Digital temperature display

Audible and visual high and low temperature alarms

2 probes; 1 in air and 1 in sample bottle

Remote alarm contacts

Adjustable operating temperature range: -15°C to -25°C

Exterior dimensions: Width: 23 <sup>5</sup>/<sub>8</sub>" , Depth: 24 <sup>1</sup>/<sub>4</sub> " , Height: 32 <sup>3</sup>/<sub>16</sub>"

Manual defrost

115V, 60 Hz, 1.2 Amps, 1/9 HP, UL/C-UL listed