Product Name:	ABT-HC-SLS-36
Description:	Standard Solid Door Laboratory Refrigerator 36 Cu. Ft. capacity
	Two swing solid doors, self-closing, right hinged, not reversible Microprocessor temperature controller
	Digital temperature display
	Audible and visual high and low temperature alarms Remote alarm contacts
	2 probes; 1 in air and 1 in sample bottle
	Adjustable operating temperature range: 1°C – 10°C
	Keyed door locks Magnetic door gasket for positive seal
	Eight adjustable shelves (adjustable in ½" increments)
	Forced draft air circulation
	Heavy duty door handles
	LED interior lights are safety shielded and switch controlled
	Swivel Casters (3 ½") Probe access port (%")
	White powder-coated steel exterior and interior
	Hydrocarbon, natural refrigerant (R290)
	One year parts and labor warranty, plus an additional four year compressor parts warranty
	Exterior dimensions: $39 \frac{5}{8}$ " W x $34 \frac{3}{4}$ " D x $81 \frac{3}{4}$ " H with casters; 78
	¹ / ₄ " H without casters
	Cycle defrost
	115V, 60 Hz, 4.5 Amps, 1/4 HP
	ETL listed
	Shipping weight: 360 lbs.