Product Name: AB7	T-HC-UCBI-0420SS
4.2 One Mic Digi Aud 2 pr Rer Adju Two Key Red Mag Lev Proi Stai Higl HFO Two Pyx Exte Mar 115 Ene UL/ ENE Ship *Th	emier Undercounter Freezer Built-In a Cu. Ft. Capacity le stainless steel swing door, right hinged* croprocessor temperature controller gital temperature display dible and visual high and low temperature alarms brobes; 1 in air and 1 in sample bottle mote alarm contacts justable operating temperature range: -15°C – -25°C for fixed shelves yed door lock cessed door handle lignetic door gasket for positive seal veling legs be access port (3/8") sinless steel exterior gh density urethane foam cabinet and door insulation C-free refrigerant (R600a), foam insulation and packaging for year parts and labor warranty xis®, Omnicell® and AcuDose RX® compatible terior dimensions: 23 ¾ "W x 24 ½" D x 33 3/8 "H anual defrost 5V, 60 Hz, 1.5 Amps, 1/8 HP ergy Consumption 0.97 kWh/day /C-UL listed liERGY STAR CERTIFIED lipping weight: 131 lbs. his product is available left hinged, product number IT-HC-UCBI-0420SS-LH