

ABS ABT-HC-UCFS-0220M Specifications and Details

1.5 Cu. Ft. Standard Undercounter Freezer Freestanding

One fixed shelf

Keyed door lock

Magnetic door gasket for positive seal

Heavy duty door handle

Leveling legs

High density urethane foam cabinet and door insulation

HFC-free refrigerant (R600a), foam insulation and packaging

One year parts and labor warranty

Mechanical thermostat

Average cycle temperature: -20°C

Exterior dimensions: Width: 18 ½", Depth: 19 ½", Height: 19 ½"

Manual defrost

115V, 60 Hz, 1/10 HP, 1.1 amps