| Product Name: | ABT-MFP-14 |
|---------------|---|
| | |
| Description: | Premier Manual Defrost Laboratory Freezer |
| | 14 Cu. Ft. capacity |
| | One solid swing door, right hinged, not reversible |
| | Four Inner "fast-freeze compartments" with doors for temperature |
| | maintenance and recovery |
| | Microprocessor Temperature Controller |
| | Adjustable operating temperature range: -15°-25°C |
| | Standard Digital Temperature Display and Alarm Module |
| | Sample simulated digital temperature display; .1° Resolution |
| | Adjustable temperature display in Celsius or Fahrenheit |
| | Audible and visual high/low temperature alarms |
| | Sensor error alarm |
| | Power failure alarm |
| | Resettable 1-point Min/Max History |
| | Remote alarm contacts |
| | Alarm validation (audible, visual, and remote alarm) |
| | Alarm mute with ring back |
| | Three fixed shelves |
| | Keyed door lock |
| | Magnetic door gasket for positive seal |
| | Heavy duty door handle |
| | Caster platform (adds 4" to height) |
| | One probe access port (3/8") |
| | White textured cabinet and door |
| | High density urethane foam cabinet and door insulation |
| | CFC-free refrigerant (R134a), foam insulation and packaging |
| | Two year parts and labor warranty; plus additional three year |
| | compressor parts warranty |
| | Exterior dimensions: 29 1/8"W x 32 3/4"D x 66"H with caster platform; |
| | 62"H without caster platform |
| | Manual defrost |
| | 115V, 60 Hz, 5 Amps, 1/3 HP |
| | UL/C-UL listed |
| | Shipping weight: 223 lbs. |