Draduat Nama	PH-ABT-HC-RFC12GA
Product Name:	FIFADI-FIO-NECIZOA
Description:	Premier Pharmacy/Vaccine Refrigerator & Freezer Combination
·	with Auto Defrost Freezer
	12 Cu. Ft. capacity (10.5 Ref./1.3 Freezer)
	Two exterior swing doors (glass door refrigerator/solid door
	freezer), right hinged, not reversible
	Microprocessor temperature controller, refrigerator and freezer
	Digital temperature display, refrigerator and freezer
	Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C25°C)
	Auto defrost cycle every 6 hours, max air temp in empty freezer:
	+4.34°C; max air temp in loaded freezer: -2.89°C. Max sample temp
	in empty freezer: -4.76°C; max sample temp in loaded freezer: -9.58°C.
	Temperature monitor device (TMD) complies with the current CDC
	guidelines, with 3 years certification of calibration, "buffered" probe
	in the product simulated solution, min/max memory. F/C switchable,
	field installable, and visual & audible temp alarm
	Audible and visual high and low temperature alarms, refrigerator
	and freezer
	2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and
	Freezer
	Remote alarm contacts, refrigerator and freezer
	Six adjustable shelves with guard rail on back to prevent product
	from being pushed against back wall, Refrigerator
	Two adjustable shelves with guard rail on back to prevent product
	from being pushed against back wall, Freezer
	Bottom footed shelf to allow air circulation, Refrigerator and Freezer
	Forced draft circulation Refrigerator Vaccine storage power cord warning label
	Keyed door locks
	Magnetic door gasket for positive seal
	Integrated door handles
	Leveling legs
	Stacking kit
	Two probe access ports (3/8"); one, refrigerator/one, freezer
	White powder coated exterior
	HFC-free refrigerant (R600a), foam insulation and packaging
	Pharmacy refrigerator/freezer toolkit and temperature logs
	Two year parts and labor warranty, plus an additional three year
	compressor parts warranty
	Exterior dimensions: 23 ¾"W x 26 ½"D x 80 ¼"H
	Auto defrost freezer / cycle defrost refrigerator
	115V, 60 Hz, 2.3 Amps refrigerator, 1 Amp freezer, 1/3 HP
	refrigerator, 1/8 HP freezer

UL/C-UL listed

Shipping weight: 316 lbs.

*Please note this unit requires 2 outlets.

DISCLAIMER: Auto defrost freezers incorporate an electric heater installed on the evaporator that assist in the removal of the frost and ice. Although the heater helps speed up the frost and ice removal, it can raise the temperatures within the freezer.