



Model **1000 + Package**

**AUTOMATIC
MANUAL**

AFFORDABLE AUDIOMETER WITH PRINTER

**Designed
For
Primary
Care
Providers**



**THE
FUTURE
IS
HERE**

Model 1000 + OTO-WITH PRINTER

Less Expensive Than Some Manual Audiometers

Reimbursable CPT Codes 92551, 92552

Electronic Medical Records Hardware Compatible

EASY TO USE – Industry's Only Five Year Warranty

Automatic Threshold Test

Ambco Audiometer Model 1000 + OtoScreen Serial Number :1197 ANSI S3/b - 1996		
AUTOMATIC HEARING TEST		DATE _____
NAME _____		ID _____
FREQUENCY	RIGHT	LEFT (EAR)
500 Hz	20	20
1000 Hz	20	20
2000 Hz	20	20
3000 Hz	25	25
4000 Hz	25	25
6000 Hz	20	30
8000 Hz	25	30
Tested by _____		

**Sample
Hard Copy
Results**

Screening Test

Ambco Audiometer Model 1000 + OtoScreen Serial Number :1197 ANSI S3/b - 1996	
NAME _____ DATE _____	
Screening (Hz) 1000 2000 3000 4000 Hz	
RIGHT EAR	Fail Pass Pass Pass
dB: 20	
1000 2000 3000 4000 Hz	
LEFT EAR	Fail Pass Pass Pass
Tested by _____	

Contact your Medical Equipment Supplier or

AMBCO 800-345-1079 www.ambco.com