

Biodex Clear-Lead™ Mobile X-Ray Barriers

Save Time – new sleek design is easy to clean!

Clear-Lead™ acrylic mobile x-ray barriers are designed for use around any imaging procedures using ionizing radiation, e.g., fluoroscopy, x-ray, CT, mammography.

The greatest advantage of these barriers over others on the market is the time they save behind the scenes. Their sleek design exposes no hardware, nooks or crannies while the opaque section of the shield is covered with stainless steel – the result, a smooth, easy-to-clean finish.

Grab-handle edges and large casters make them effortless to maneuver. Biodex Clear-Lead™ Mobile Barriers offer large, durable and shatter resistant protection wherever it's needed.

Features

- Meets all infection control criteria
- Sleek design is easy to clean
- Durable and shatter resistant protection
- Effortless maneuverability
- 0.5 mm lead equivalency window
- 0.8 mm lead stainless steel panel
- Superior stability

Specifications

042-580 Clear-Lead Mobile Barrier, Standard Window

- Dimensions:
 - Overall: 33" w x 73.62" h (84 x 187 cm)
 - Window: 30" w x 22" h (76 x 56 cm)
 - Opaque Panel: 30" w x 46.25" h (76 x 117 cm)
 - Shielded Area: 30" w x 68.25" h (76 x 173 cm)
- Shielding:
 - Window: 0.5 mm lead equivalency
 - Opaque Panel: 0.8 mm lead
- Leg Depth: 10.5" (26.7 cm)
- Casters: Four hospital grade; two locking, two non-locking
- Weight: 103 lb (47 kg)
- Warranty: One year parts and labor

042-582 Clear-Lead Mobile Barrier, Full Window

- Dimensions: overall: 33" w x 73.62" h (84 x 187 cm)
 - Window: 30" w x 46" h (76 x 116.8 cm)
 - Opaque Panel: 30" w x 22.25" h (76 x 56.5 cm)
 - Shielding Area: 30" w x 68.25" h (76 x 173 cm)
- Shielding:
 - Window: 0.5 mm lead equivalency
 - Opaque Panel: 0.8 mm lead
- Leg Depth: 10.5" (26.7 cm)
- Casters: Four hospital grade; two locking, two non-locking
- Weight: 97 lbs. (44 kg)
- Biodex Warranty: One year parts and labor