



ClearSpec

clearspecmedical.com



Empowering Better Visualization



ClearSpec® is the New Standard of Care

- Better visualization for every exam
- Improves physician workflow
- Maximizes patient comfort

New ClearSpec® LED Illuminator*

- Reposable high power LED light source
- Fits into the handle of the ClearSpec® speculum
- Provides more than 2 hours of light
- No more cords or recharging batteries




*Patent pending

Advanced Speculum Design

ClearSpec® is scientifically proven by an IRB randomized clinical study to improve visualization for every exam¹


The Problem

Standard Speculums only provide single-plane (superior/inferior) tissue retraction



The ClearSpec® Solution

ClearSpec® Speculums offer multi-planar tissue retraction




Improved visualization with side wall retraction and high powered LED illumination

The ClearSpec® Advantage



View of the cervix using a standard speculum



View of the cervix using ClearSpec®

MFG#	Item Descriptions	Additional Features	UOM
102	ClearSpec® Single Use Vaginal Speculum	Small, Graves Style, Non-sterile	25 ea/bx, 4 bx/cs
103	ClearSpec® Single Use Vaginal Speculum	Medium, Graves Style, Non-sterile	25 ea/bx, 4 bx/cs
104	ClearSpec® Single Use Vaginal Speculum	Large, Graves Style, Non-sterile	25 ea/bx, 4 bx/cs
105	ClearSpec® LED Speculum Illuminator	Battery Operated Light Source	1 ea/bx, 12 bx/cs

Not made with natural rubber latex.