

CONSTRUCTION –


- High impact “no break” polyethylene
- One-piece roto-molded construction
- Manufactured to ISO Quality Standards

FEATURES –

- Limited Lifetime Warranty
- Bedside has integrated rails to help contain items on top surface
- Bed and cabinet sides fit flush to walls
- Bed features surface grooves for fluid management
- Single side bed storage can be positioned left or right facing
- Stain-resistant, seamless and easy to clean
- Corners and edges rounded for increased safety
- 100% recyclable material

OPTIONS –

- Enclosed bottom
- Bolt-to-floor brackets
- Patient restraint module



	FLBD3680 Bed	FLBDS3680 Bed, side storage	FL100 Bedside Cabinet	FRCG702 Two Shelf Wall Mount Unit	FRCG703 Three Shelf Wall Mount Unit
Width	39.5	39.5	20	29	29
Depth	83.5	83.5	18	14	14
Height	20	20	23.4	22	51
Deck Height	18	18	-	-	-
Top Height	-	-	22.4	-	-
Storage Width	-	6.5/20	15	-	-
Storage Depth	-	16	15.75	-	-
Storage Height	-	11	5/8.5	-	-
Cubic Feet	33	33	5	5	12
Ship Weight	146 lbs	155 lbs	26 lbs	23 lbs	49 lbs



	SF36805 5" Mattress	SF36805-PW 5" mattress, built in pillow
Width	80	80
Depth	36	36
Height	5	5
Cubic Feet	14	14
Ship Weight	18 lbs	18 lbs