

FR566G

LUMEX[®]
by graham-field

Deluxe Clinical Care Recliner is ideal for clinical and acute care settings. It provides the best value while being engineered for maximum patient comfort, durability, and ease of use. It features a contoured back to support the patient comfortably in all three positions: upright, first recline and second recline. Red color-coded side pedals quickly put the patient into Trendelenburg position while allowing the caregiver to maintain eye contact with the patient. Fold down side tables with smooth surfaces and removable side panels for easy cleaning.



CLINICAL
CARE



**400
POUND**¹
WEIGHT CAPACITY
EVENLY DISTRIBUTED

PACU	Rehab	Oncology	IV Therapy	Exam Room	Patient Room	Emergency Room	Dialysis (Facility)	Dialysis (In Home)	Nocturnal Dialysis	Physical Therapy	Pain Management	Post Operative Care	Same Day Surgery / ASC
■	■	■	■	■	■	■	■	■	■	■	■	■	■

AREAS OF USE

Options:

* Model	Wallsaver Armrest Caps	Fully Upholstered Arms	Heat and Massage	Heat Only
FR566Gxxx	■			
FR566GHxxx	■		■	
FR566GHOxxx	■			■
FR566GxxxPA		■	■	
FR566GHxxxPA		■	■	
FR566GHOxxxPA		■		■

† xxx REPRESENTS FABRIC COLOR

 **FR566G** Product Info



Positions: THREE POSITION + CAREGIVER-ACTIVATED TRENDELENBURG

Three Positions for Patient Comfort: Upright, First Recline and Second Recline
Fourth Position: Trendelenburg (or Shock) Position, Caregiver-Activated Using Side Pedal



Features:

- Wide and contoured seat and back provides an enhanced level of patient comfort
- Higher seat & back height allows caregiver to more comfortably assist patient
- Integrated leg and footrest for exceptional comfort
- 5 inch Tente® Swivel casters for ease of maneuverability: 3 Total Locking and 1 Directional Locking
- Wide and soft armrest for patient comfort
- Wall saver armrest caps increase durability and protect walls from marring; – also available with fully upholstered arms
- Two self-storing, fold down, side tables with smooth surface for easy cleaning
- Smooth and removable side panels for easy cleaning and maintenance
- Rear push handle for improved patient transport



FOLD DOWN, FLUSH-MOUNT SIDE TABLES



SOFT, WIDE ARMREST




TENTE® CASTERS



REMOVABLE SIDE PANELS



REAR PUSH HANDLE



INTEGRATED LEG AND FOOTREST

Accessories:

★ sold separately

	item no.†
Headrest Pillow	FR5659xxxUS
Headrest Skirt	HRSxxx
Headrest Cover	HRC566xxx
4 Inch Extension Legs ⁶	FR5699G
Swing-away Folding Table	SAFT10
IV Pole	2101A
IV Mount	5782RG

Specifications:

Overall Width		Overall Width side tables up		Overall Depth		Overall Depth trendelenburg		Overall Height		Arm Width		Width between arms		Arm Height from floor		Seat Width		Seat Depth		Seat Height		Back Width		Back Height above seat		Back Angle upright legrest extended fully reclined	Weight ¹ Capacity EVENLY DISTRIBUTED			
inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†		lb	kg†		
31.25	79.38	51.50	130.81	37.50	95.25	75.00	190.50	49.50	125.73	4.50	11.43	21.50	54.61	28.50	72.39	20.50	52.07	19.00	48.26	22.25	59.52	27.50	69.85	32.00	81.28	120	130	150	400	181.44

* Ships Freight: request additional shipping quote † xxx: in the product number is the placeholder for fabric choice ¹ maximum weight capacity - evenly distributed - refer to user manual
⁶ adjust height specs: +4 inches (+10.16 cm) † standard rounding method used ★ options available at additional cost - accessories sold separately

Available Colors:

COLOR OPTIONS AND FABRIC SWATCHES ON WEBPAGE  **HC Fabric Choices**



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.
Tente is a registered trademark of TENTE International GmbH.

graham field

GF09000.08Rev125