

Deluxe Tilt Overbed Table

- Easy assembly.
- Chrome-plated and welded tubular steel column and base make the table durable and stable.
- 1 1/2" hooded ball-bearing casters provide easy maneuverability across all floor surfaces.
- Table top adjusts 13" in infinite positions and is constructed of wood-grained oak Melamine laminate on particle-board core.
- Raised, T-molded plastic edge.
- Chrome-plated "H" base and column for aesthetics and durability.
- Designed to meet the requirements of HCPCS code: E0274.
- Supports up to 30 lb of dispersed weight.

GF8905-1A 1 ea

ITEM	Table Height	Table Width	Table Depth	Base Width	Base Depth	Weight	Dispersed Weight Capacity
GF8905-1A	30" – 43.8"	30"	15.2"	26.2"	16.5"	21 lb	30 lb



GF8905-1A



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com