



**Advantages**

- EMR Connectivity**  
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter<sup>®</sup> Professional Connectivity Solutions
- Measuring Tape Included**  
Easy to read measuring tape allows user to measure both weight and height simultaneously
- Display Head 17" Above Scale Base**  
Display head is positioned 17" above the scale base to allow for an unobstructed view of the weight measurement
- Value-Added Features**  
Operating functions (LB/KG Conversion, LB/KG Lock Out, Zero Out/Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost

**Health o meter®**  
P r o f e s s i o n a l



Hygienic,  
easy clean tray

Measuring Tape  
Included

Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
522KL	26.25" x 17.5" x 14.25" / 3.79	24.25 lb	892076002541



**Specifications**

- Capacity: 50 lb / 23 kg
- Resolution: 0 - 20 lb 0.2 oz / 20 - 50 lb / 0.5 oz; 0 - 9 kg / 5g; 9 - 23 kg / 10 g
- Connectivity via USB
- 1" LCD display
- Functions: LB / KG Conversion, LB /KG Lock Out, Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- Tray Size (w x d x h): 24 1/8" x 14 1/2" x 2 5/8" / 613 mm x 368 mm x 67 mm
- Measuring Tape: 0" - 23" / 0 cm - 58 cm
- 6 AA batteries included, Optional 120V adapter (order # ADPT31)
- 2 Year Limited Warranty