



Advantages

- ✓ **EMR Connectivity**
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- ✓ **Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, Zero Out/Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- ✓ **1" LCD Display**
The large display shows the user an easy-to-read weight measurement
- ✓ **Convenient Weighing Seat**
Plastic, easy-to-clean seat limits infant movement while accurately calculating weight measurement



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
524KL	26.375" x 17.25" x 20.625" / 5.43	26.5 lb	892076002541



Specifications

- Capacity: 50 lb / 23 kg
- Resolution: 0 - 20 lb 0.2 oz / 20 - 50 lb 0.5 oz 0 - 9 kg / 5g; 9 - 23 kg / 10 g
- Connectivity via USB
- 1" LCD display
- Functions: LB / KG Conversion, LB /KG lock out, Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- Seat Size: (w x d x h) 24 1/4" x 14 1/4" x 18 1/2" / 616 mm x 362 mm x 470 mm
- 6 AA batteries included, Optional 120V adapter (order # ADPT31)
- 2 Year Limited Warranty