

2500KL
Digital Wheelchair Ramp



2500CKL
Digital Wheelchair Ramp Scale
with Fold Away Seat



Advantages

- ✓ **EMR Connectivity**
Scale interfaces with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- ✓ **Easy Portability**
The scale has two wheels that provide easy mobility and allow for convenient placement
- ✓ **Swivel Display Head**
For the convenience of the administrator, the display head rotates and pivots allowing line-of-sight viewing and information privacy
- ✓ **On-Screen Help Menu**
Guides users through operating functions and provides easy access to multiple settings
- ✓ **Everlock™ (patent pending)**
Option that allows users to permanently lock the scale in either LB or KG to improve safety.
- ✓ **Interactive Display**
High-contrast color screen provides a user-friendly interface
- ✓ **Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option) included at no additional cost



2500CKL / 2500CKG



Specifications

- Capacity: 1,000 lb / 454 kg
- Seat Capacity: 419 lb / 190 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size (w x h x d):
29 3/8" x 28 1/4" x 2 3/8" / 746 mm (w) x 718 mm (d) x 60 mm (h)
- Seat Size (for 2500CKL, 2500CKG or CHAIRACC):
18" (w) x 15" (d) x 18 1/2" (h) / 457 mm (w) x 381 mm (d) x 470 mm (h)
- Ramp Size: 28 1/4" (w) x 9 7/8" (d) / 714 mm (w) x 248 mm (d)
- Product Footprint (w x h x d):
40 1/4" x 34 1/2" x 45 1/2" / 1022 mm / 876 mm / 1156 mm
- Product Weight:
- Display: 1 1/2" / 38 mm High-Contrast Color TFT-LCD Screen
- Wheels: 2
- Power Source: 100-240V adapter included, 6 D batteries (not included)
- Connectivity: USB, optional Bluetooth®
- Functions: LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option, Everlock™
- 2 year limited warranty

Model	Description
2500KL	Digital wheelchair ramp scale, pounds and kilograms
2500KG	Digital wheelchair ramp scale, kilograms only
2500CKL	Digital wheelchair ramp scale with fold away seat, lb & kg
2500CKG	Digital wheelchair ramp scale with fold away seat, kilograms only
SS-2500KL	ScaleSurance 2 year extended warranty
SS-2500CKL	ScaleSurance 2 year extended warranty
CHAIRACC	Optional fold away chair attachment for 2500KL
ADPT30	Adapter Included, item # for replacement or spare
B2500RAMP	Optional secondary ramp for 2500KL & 2500CKL
BLUETOOTHONE-P	Wireless thermal printer