

Vascular and Obstetric Probes

APPG probe



4 MHz Vascular Probe



5 MHz Vascular Probe



8 MHz Vascular Probe



10 MHz Vascular Probe



EZ8 8 MHz Vascular Wide Beam Probe



2 MHz Obstetric Probe



3 MHz Obstetric Probe




Warranty

1 year for Probe

Configurations

PPGA1	Arterial Plethysmography (APPG) sensor for the toe
VP4XS-USA	Detection of deep lying vessels
VP5XS-USA	Edematous limbs and deep lying vessels
VP8XS-USA	Detection of peripheral vessels and calcified arteries
VP10XS-USA	Detection of smaller vessels in specialist superficial applications
EZ8XS-USA	Easy location of vessels
OP2XS-USA	Fetal Heart Rate detection
OP3XS-USA	Fetal Heart Rate detection in early gestation



Contact:

Arjo, Inc. 2349 W. Lake St., Suite 250, Addison, IL 60101
Phone: +1 (800) 323 1245 Fax: +1 888 594 2756
Web: www.huntleigh-healthcare.us Email: Diagnosics.USADD@arjo.com

HUNTLEIGH

A MEMBER OF THE ARJO FAMILY