

## Precision Medical Easy Go Vac Aspirator

Power & long life on the go in a battery operated suction unit!

<b>Specifications</b>	
Dimensions (in Carry Bag)	Length 12" (30.3 cm), Width 9 1/2" (24.1 cm), Height 10 3/16" (25.9 cm)
Weight (with Carry Bag and all Accessories)	9.6 lbs. (4.4 kg)
Shipping Weight	11 lbs. (5 kg)
Electrical	Aspirator unit: 12 Volts DC Battery Charger/Power Supply: 100-240 Volts AC
Operating Temperature Range	50°F to 104°F (10°C to 40°C)
Storage Conditions	Temperature: -4°F to 104°F (-20°C to 40°C) Humidity: Max 95% non-condensing
Vacuum Range	50 mmHg to 535 mmHg (2 inHg to 21 inHg)
Flow	14 lpm (unrestricted flow, without suction canister, running on AC power)
Power Consumption	20 watts
Vacuum Buildup	0 to 300 mmHg in 3.1 sec. (with 800cc Suction Canister and AC Power Cord connected) 0 to 300 mmHg in 3.5 sec. (with 800cc Suction Canister connected and running on Battery)
Suction Canister	800cc Disposable with Float Shut-off