


Product Comparison Chart

Features	641-002	641-003	641-004	641-005	630-003	630-004	630-005	630-006	230-001	230-002	230-003
Patient weight capacity	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.	450 lbs.
Minimum seat height	19"	19"	22.5"	22.5"	19"	19"	22.5"	22.5"	22.5"	22.5"	26"
Maximum seat height	40"	40"	43.5"	43.5"	40"	40"	43.5"	43.5"	40"	40"	43.5"
Headrest upholstery dimension	9"	9"	9"	9"	12"	12"	12"	12"	18"	18"	18"
Back upholstery dimension	21.5"	21.5"	21.5"	21.5"	24" - 26"	24" - 26"	24" - 26"	24" - 26"	27"	27"	27"
Seat upholstery dimension	23.5"	23.5"	23.5"	23.5"	28"	28"	28"	28"	27"	27"	27"
Foot upholstery dimension	21.75"	21.75"	21.75"	21.75"	18"	18"	18"	18"	18"	18"	18"
Foot control	•	•	•	•	•	•	•	•	• opt	• opt	• opt
Hand control	•	•	•	•	•	•	•	•	•	•	•
QuickChair™/Auto-return	•	•	•	•	•	•	•	•			
Crash-avoidance safety system	•	•	•	•	•	•	•	•			
Premium Comfort System™	•	•	•	•	•	•	•	•			
Hospital grade electrical outlets		•			•	•				•	
Side rail adjustment	•	•	•	•							
Flip-out foot extension	•	•	•	•							
Stirrups/Stainless steel treatment pan					•	•	•	•	•	•	•
Knee-chest positioning									•	•	•
Programmability		•		•		•		•			
Rotation			•	•			•	•			•
Paper roll holder	• opt	• opt	• opt	• opt	•	•	•	•	•	•	•
Wireless controls	• opt	• opt	• opt	• opt	• opt	• opt	• opt	• opt			
Ultra-Premium upholstery	• opt	• opt	• opt	• opt	• opt	• opt	• opt	• opt			
Upholstery colors	27	27	27	27	27	27	27	27	8	8	8

opt = optional


Midmark Premium Solid Upholstery Colors



Midmark Premium Pattern Upholstery Colors



Midmark Ultra-Premium Solid Upholstery Colors



Ritter Solid Upholstery Colors



For product inquiries, please contact Midmark Medical Customer Service at 1-800-MIDMARK.
 If you have questions regarding the promotion, please contact the Promotions Department, ext. 7575.
 For more information or a demonstration, contact your Midmark or Ritter Dealer, or Midmark, 60 Vista Drive, P.O. Box 286, Versailles, Ohio 45380-0286. Call: 1-800-MIDMARK Fax: 1-800-365-8631
midmark.com

