

Seiler Instrument Microscope Division

Products:

SeilerScope

- Reversed four-hole nose piece
- Inclined binocular tube (f=1x)
(Seidentoph type)
- 2 - 10x widefield eyepieces (18)
- Abbe condensor 1.25 NA with iris
(Centerable)
- Halogen illumination (6V/20W)
(Replacement Bulb #64250)
- Achromat objectives (4x, 10x, 40x,
and 100x Oil)
- Graduated, right hand mechanical
stage
- Sturdy base with built-in transformer
- Focusing tension adjustment knob
- Rack and pinion mount
- Dust cover



The SeilerScope

The newest member of the Seiler family of microscopes offers features and benefits of higher priced microscopes at a much lower cost. A real bargain, ideal for the office.