



Specialists



Cardiology



Vascular Surgery



Interventional Radiology



Rapid ABI-Q



Single-Level



Multi-Level



TBI w/PPG

CLINIC PLACEMENT

EXAMS PERFORMED



American-Made Quality



Includes labor, parts, and accidental breakage



Automatic Cuff Selector



Push-Button Remote



simpleABI Roll Stand



ABI-Q Rapid Screen



DICOM Compatible

simpleABI


PAD testing has never been easier thanks to our complete family of user friendly simpleABI systems that provide the information and reports you need quickly and intuitively. All of our vascular ABI systems are Made in the USA and have an industry-leading 2 year warranty.

ABI-500CL Exams

- ABI-Q Rapid Screen
- Single-Level
- Multi-Level
- TBI with PPG

ABI-500CL Includes:

- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - 2 x 12cm EZ-Clean Vasocuffs
 - 2 x 2.5cm disposable digit cuffs
- Two-Year Warranty
 - Includes all parts, labor, and accidental breakage



American-Made Quality

- All ABI and Doppler systems are manufactured outside of Denver, Colorado using the highest-quality components available
- Our systems are built to last and endure years of use in the busiest of clinical settings
- Our industry-leading 2-year warranty covers all labor, parts, and accidental breakage

Newman Medical History & Commitment

- The Newman family has been involved with ultrasound diagnostics since the 1970s
- Newman Medical is committed to ongoing product development and bringing high-quality vascular diagnostic products to market
- Newman Medical has a relentless focus on service and establishing long term relationships with our customers

First-Class Support

- Newman Medical boasts the industry's best and most tenured customer service team with 25+ years' average experience
- Newman's extensively trained field representatives are always available to support physician needs
- Newman Medical offers both on-line and in-person clinical support to physicians


ABI-500CL




Exams Performed



ABI-Q Rapid Screen



Single-Level



Multi-Level



TBI with PPG

Ordering Information

System

SKU
ABI-500CL

DESCRIPTION
simpleABI Multi-Level Cuff-Link System

Accessories

SKU
ABI-120
STND-300
ACC-190
D5
SFTW-100
ABI-WARR

DESCRIPTION
Microsoft® Windows® PC for use with simpleABI Systems
simpleABI Roll Stand for Cuff-Link Systems
Large simpleABI carry bag
5Mhz vascular probe only
DICOM software integration
1 year additional warranty for simpleABI systems

Replacement Cuffs

SKU
CUFF-200
CUFF-210
CUFF-210L
CUFF-230
CUFF-250
CUFF-220D
CUFF-240D
CUFF-SET-242

DESCRIPTION
EZ Clean Vascular Cuff, 10cm (ankle, small arm)
EZ Clean Vascular Cuff, 12cm (arm, ankle, calf)
EZ Clean Vascular Cuff, 12cm long (above knee)
EZ-Clean Vascular Cuff, 17cm Contoured (thigh)
EZ-Clean Vascular Cuff, 7cm long (metatarsal)
Digit Cuff, 2.5cm - Latex-Free
Digit Cuff, 1.9cm - Latex-Free
EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk

ABI-500CL Technical Specifications

Height: 107 cm (42 inches)

Weight: 5.5 kgs (12 pounds)

Level of Protection against electrical shock: Type B Applied Part, Class II Equipment

Comply with the following Standards: IEC60101-1, IEC60601-2, IEC60601-2-37

Operating Temperature: 10° ~ 40°C (50° ~ 104°F)

Operating Humidity: 30% ~ 75%

Transport/Storage Temperature: -20° ~ 50°C (-4° ~ 122°F)

Transport/Storage Humidity: 5% ~ 90%, non-condensing

Battery Life: 10 hours

GTIN: 00852785008616

REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

