




Primary Care



Podiatry



Wound Care



Single-Level



TBI w/PPG

CLINIC PLACEMENT

EXAMS PERFORMED



American-Made Quality



Includes labor, parts, and accidental breakage



DICOM Compatible

simpleABI

PAD testing has never been easier thanks to our complete family of user friendly simpleABI systems that provide the information and reports you need quickly and intuitively. All of our vascular ABI systems are Made in the USA and have an industry-leading 2 year warranty.

ABI-300 Exams

- Single-Level
- TBI with PPG

ABI-300 Includes:

- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - 2 x 12cm EZ-Clean Vasocuffs
 - 2 x 2.5cm disposable digit cuffs
- Two-Year Warranty
 - Includes all parts, labor, and accidental breakage



American-Made Quality

- All ABI and Doppler systems are manufactured outside of Denver, Colorado using the highest-quality components available
- Our systems are built to last and endure years of use in the busiest of clinical settings
- Our industry-leading 2-year warranty covers all labor, parts, and accidental breakage

Newman Medical History & Commitment

- The Newman family has been involved with ultrasound diagnostics since the 1970s
- Newman Medical is committed to ongoing product development and bringing high-quality vascular diagnostic products to market
- Newman Medical has a relentless focus on service and establishing long term relationships with our customers

First-Class Support

- Newman Medical boasts the industry's best and most tenured customer service team with 25+ years' average experience
- Newman's extensively trained field representatives are always available to support physician needs
- Newman Medical offers both on-line and in-person clinical support to physicians

ABI-300



Exams Performed



Single-Level



TBI with PPG

Ordering Information

System

SKU	DESCRIPTION
ABI-300	simpleABI Single-Level Cuff-Link System with TBI

Accessories

SKU	DESCRIPTION
ABI-120	Microsoft® Windows® PC for use with simpleABI Systems
STND-130	simpleABI Roll Stand for Manual Systems
ACC-175	simpleABI-300 carry bag
PUMP-100	Manual cuff inflator

Replacement Cuffs

SKU	DESCRIPTION
CUFF-200	EZ Clean Vascular Cuff, 10cm (ankle, small arm)
CUFF-210	EZ Clean Vascular Cuff, 12cm (arm, ankle, calf)
CUFF-210L	EZ Clean Vascular Cuff, 12cm long (above knee)
CUFF-220D	Digit Cuff, 2.5cm - Latex-Free
CUFF-240D	Digit Cuff, 1.9cm - Latex-Free
CUFF-SET-242	EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk

ABI-300 Technical Specifications

Weight:	3.6 kgs (8 pounds)
Level of Protection against electrical shock:	Type B Applied Part, Class II Equipment
Comply with the following Standards:	IEC60101-1, IEC60601-2, IEC60601-2-37
Operating Temperature:	10° ~ 40°C (50° ~ 104°F)
Operating Humidity:	30% ~ 75%
Transport/Storage Temperature:	-20° ~ 50°C (-4° ~ 122°F)
Transport/Storage Humidity:	5% ~ 90%, non-condensing
Battery Life:	800, 1 minute exams
GTIN:	00852785008647

REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

