

UNICO G380 Infinity Microscopes

Model Configuration

MODEL	HEAD	EYEPIECE	OBJECTIVE
G386 G386PL	Binocular	High-Eyepoint WF10X/18	Infinity Achromatic 4X/10X/40X/100X(oil) Infinity Plan Achromatic 4X/10X/40X/100X(oil)
G388 G388PL	Dual Binocular	High-Eyepoint WF10X/18	Infinity Achromatic 4X/10X/40X/100X(oil) Infinity Plan Achromatic 4X/10X/40X/100X(oil)

	G386	G386PL	G388	G388PL
Head	Binocular	Binocular	Dual-Binocular	Dual-Binocular
Eyepiece	10X WF(2)Hi-EP	10X WF(2)Hi-EP	10X WF(4)Hi-EP	10X WF(4)Hi-EP
Objectives	4,10,40,100X Infinity Achromat	4,10,40,100X Infinity Plan	4,10,40,100X Infinity Achromat	4,10,40,100X Infinity Plan
Condenser	NA 1.25	NA 1.25	NA 1.25	NA 1.25
Diaphragm	Iris	Iris	Iris	Iris
Stage	Mechanical	Mechanical	Mechanical	Mechanical
Illumination	Halogen	Halogen	Halogen	Halogen
Bulb	6 Volt/20 Watt	6 Volt/20 Watt	6 Volt/20 Watt	6 Volt/20 Watt
Focusing	Coaxial	Coaxial	Coaxial	Coaxial