

UNICO G380 Series Microscopes

Model Configuration

MODEL	HEAD	EYEPIECE	OBJECTIVE
G380 G381 G382 G383	Binocular Monocular Dual View Trinocular	High-Eyepoint WF10X/18	Achromatic 4X/10X/40X/100X(oil)
G380-SP G381-SP G382-SP G383-SP	Binocular Monocular Dual View Trinocular	High-Eyepoint WF10X/18	Semi-Plan 4X/10X/40X/100X(oil)
G380-PL G381-PL G382-PL G383-PL	Binocular Monocular Dual View Trinocular	High-Eyepoint WF10X/18	Plan Achromatic 4X/10X/40X/100X(oil)

	G380	G380SP	G380PL	G381	G381SP	G381PL
Head	Binocular	Binocular	Binocular	Monocular	Monocular	Monocular
Eyepiece	10X WF(2)Hi-EP	10X WF(2)Hi-EP	10X WF(2)Hi-EP	10X WF Hi-EP	10X WF Hi-EP	10X WF Hi-EP
Objectives	4,10,40,100X Achromat	4,10,40,100X Semi-Plan	4,10,40,100X Plan	4,10,40,100X Achromat	4,10,40,100X Semi-Plan	4,10,40,100X Plan
Condenser	NA 1.25	NA 1.25	NA 1.25	NA 1.25	NA 1.25	NA 1.25
Diaphragm	Iris	Iris	Iris	Iris	Iris	Iris
Stage	Mechanical	Mechanical	Mechanical	Mechanical	Mechanical	Mechanical
Illumination	Halogen	Halogen	Halogen	Halogen	Halogen	Halogen
Bulb	6 Volt/20 Watt	6 Volt/20 Watt	6 Volt/20 Watt	6 Volt/20 Watt	6 Volt/20 Watt	6 Volt/20 Watt
Focusing	Coaxial	Coaxial	Coaxial	Coaxial	Coaxial	Coaxial

	G382	G382SP	G382PL	G383	G383SP	G383PL
Head	Dual Head	Dual Head	Dual Head	Trinocular	Trinocular	Trinocular
Eyepiece	10X WF(2)Hi-EP	10X WF(2)Hi-EP	10X WF(2)Hi-EP	10X WF(2)Hi-EP	10X WF(2)Hi-EP	10X WF(2)Hi-EP
Objectives	4,10,40,100X Achromat	4,10,40,100X Semi-Plan	4,10,40,100X Plan	4,10,40,100X Achromat	4,10,40,100X Semi-Plan	4,10,40,100X Plan
Condenser	NA 1.25	NA 1.25	NA 1.25	NA 1.25	NA 1.25	NA 1.25
Diaphragm	Iris	Iris	Iris	Iris	Iris	Iris
Stage	Mechanical	Mechanical	Mechanical	Mechanical	Mechanical	Mechanical
Illumination	Halogen	Halogen	Halogen	Halogen	Halogen	Halogen
Bulb	6 Volt/20 Watt	6 Volt/20 Watt	6 Volt/20 Watt	6 Volt/20 Watt	6 Volt/20 Watt	6 Volt/20 Watt
Focusing	Coaxial	Coaxial	Coaxial	Coaxial	Coaxial	Coaxial