



DIAGNOSTIC CARDIOLOGY SUITE SPIROMETRY

Make spirometry simple with a single point of connectivity and a common interface between spirometry and ECG modalities. And with one device capable of two exams, the Diagnostic Cardiology Suite allows you to deploy, maintain and connect one secure solution for simple exams.



Simple and Secure: A solution that fits into your existing workflow.

SIMPLE

Send spirometry data to your EMR to benefit from a common workflow, user interface and setup.

SECURE

The Diagnostic Cardiology Suite Spirometry uses the same security features as the ECG modality to control user access, assign user roles and view audit trails.

ACCURATE

An incentive screen and other exam features can help you instantly verify the quality and variability of each spirometry effort. Repeatability criteria provides feedback and offers opportunities for real-time exam coaching.



SPIROMETER

Our spirometer uses disposable, single-use transducers to help reduce the risk of cross contamination between exams. Integrated ambient sensors enter temperature, pressure or humidity values before an exam, so you don't have to.