



Medical with HOTV: For the complete vision screening needs preschool through adult. Ideal for School, Athletic Physicals, Employment Physicals, and General Check-Ups.



Slide #1 (2000-026) Far Point
Acuity Binocular HOTV Test (20/70)



Slide #7 (2000-019) Far Point
Color Perception (Tumbling "E")




Slide #2 (2000-030) Far Point
Acuity Right Eye & Left Eye Monocular
HOTV Test (20/30)



Slide #8 (2000-022) Near Point
Muscle Balance
(Combination Lateral & Vertical Phoria)



Slide #3 (2000-031) Far Point
Acuity Right Eye & Left Eye Monocular
HOTV Test (20/30)




Slide #9 (2000-003) Far Point
Acuity Right Eye & Left Eye Monocular
or Both Eyes Simultaneously (20/200-20/20)



Slide #4 (2000-032) Far Point
Acuity Right Eye & Left Eye Monocular
HOTV Test (20/30)




Slide #10 (2000-012) Far Point
Lateral Phoria (1 Diopter Increments)




Slide #5 (2000-004) Far Point
Acuity Right Eye & Left Eye Monocular
Right Eye & Left Eye Plus Lens Test (20/200-20/20)



Slide #11 (2000-025) Far Point
Vertical Phoria (1/2 Diopter Increments)



Slide #6 (2000-020) Far Point
Muscle Balance
(Combination Lateral & Vertical Phoria)



Slide #12 (2000-007) Near Point
Acuity Right Eye & Left Eye Monocular
or Both Eyes Simultaneously (20/100-20/20)